

COTL1-coactosin-like 1 (Dictyostelium) Gene

Reference: PTXBC042970

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CLP; coactosin-like protein; coactosin-like 1; coactosin like F-actin binding protein 1		
For research use only.			

GENE DESCRIPTION	coactosin-like 1 (Dictyostelium)
GENE ID	23406
ORIGIN SPECIES	Human
ACCESSIONS	BC042970
VECTOR	pDONR223
NCBI SYMBOL	COTL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaccaagatcgacaagaggcttgcggggcggtacaaacctgtgcgcgacgacggctcggcgctcatctgggtgacttttaaatatgacggctccaccatcggtccacggcgagcaggagcgggtaccagcacttcaccagcagtcacagatgacgtcaggtgtttgc
cttcgtgccttcaccaccggggatgccatgagcaagaggtccaagtttgcctcatcagctggatcgggtgagaacgtcagcgggctgcagcgcgccccaaacgggacggacaagaccctggtgaaggagggtcgtacagaatttcgtaaggagtttgtatcagtgatcgaaggag
ctggaggaagatttcacagagcagctgaagaaggcggggggagccaattacgacgccagacggagtaa