

## UBL4A-ubiquitin-like 4A Gene

Reference: PTXBC043346

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DX254E; DXS254E; G6PD; GDX; GET5; MDY2; TMA24; UBL4; ubiquitin-like protein 4A; ubiquitin-like 4; ubiquitin-like protein GDX; ubiquitin like 4A		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	ubiquitin-like 4A
<b>GENE ID</b>	8266
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC043346
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	UBL4A
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

atgcagctgacggtgaaggcctgcaggccgcgagtcagcctgcaggtgccagaggacgagctgtgtccacgctgaagcagctggtctccgagaagctgaacgtcccagtcgccagcagcggctgctgtcaaggccaagccctggcagatgggaaacgactctcgatt  
atagcctcgggcccactccaagctcaacctagtggtcaaacctggagaaggtgctactagaagaagcgaggcccagagctggccgactccccacccccgaggtctgagctgatctccaagctctggcccgccactcagtcggcagatgccagcagggctctggaac  
agctcagagggattacgagaggtccctgagtcctgacgctggagcagatcgaacgggtggccagccctctctgaccctgaagtgactgagacaatggagaagggctctccaatag