

PTCH1-patched homolog 1 (Drosophila) Gene

Reference: PTXBC043542

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BCNS; HPE7; NBCCS; PTC; PTC1; PTCH; PTCH11; protein patched homolog 1; PTCH protein +12b; PTCH protein +4'; PTCH protein -10; PTCH protein -3,4,5; patched 1		

For research use only.

GENE DESCRIPTION	patched homolog 1 (Drosophila)
GENE ID	5727
ORIGIN SPECIES	Human
ACCESSIONS	BC043542
VECTOR	pDONR223
NCBI SYMBOL	PTCH1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggaaggctactggccgaaagcgcgctgtggctgagagcgaagttcagagactcttattaaactgggtgttacatcaaaaaactgcggcaagttctgtgtggcctctcatattggggccttcgcggtgggattaaaagcagcgaacctcgagaccaacgtggagga
gctgtgggtggaagttggaggacgagtaagtcgtgaattaaattatactgccagaagattggagaagaggtatgttaatcctcaactcatgatacagaccctaaagaagaaggtgctaagtctgaccaagaagcgtcctacaacactggactggcactccaggccagccgt
gtccatgtatacatgtacaacaggcagtggaattggaacattgtgttacaatcaggagagcttatcacagaacaggttacatggatcagataatagaatactttaccctgtttgattattacacctttggactgcttctgggaagggcgaaattacagtctgggacagcatatcctctgta
a