

MGC50273-MGC50273 protein Gene

Reference: PTXBC043584

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			

For research use only.

GENE DESCRIPTION	MGC50273 protein
GENE ID	408029
ORIGIN SPECIES	Human
ACCESSIONS	BC043584
VECTOR	pDONR223
NCBI SYMBOL	MGC50273
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttcagcctgagagtgagagcagggccagagacactagctctgctcaggtgcagaatacagcgalttcagtgccagctgttgagccagtgccagcaccagggcagattcccctccagggacagcgtggagctagaggaagctccagagccctctgcccgtgcc
ctgggactgcccagaccagcccagtgaggagctgaccttcaatggcacctctctagagccaccggcctcagccctggagctgaaagtgtggctggagctagaggtggcagagaggggtggccagcacagctccagccagcagctcccactgctcccagctgggacag
tggaagctatggagcagagaccaggtttgcaatcggctcctctgctcactggagagggacttctcattcagcagagcagcagccctgctgctgaagggcctgctgactgctgctggggctgtttgctgctcctcagagaggtgctggagagcaaaaaaggagcctgtgagc
aggggtccagcaggtcctctgctcccagggcagctcctccaggcaggttggcctcagctgggcatctggccalttggccctaa