

PRND-prion protein 2 (dublet) Gene

Reference: PTXBC043644

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DOPPEL; DPL; PrPLP; dJ1068H6.4; prion-like protein doppel; prion gene complex, downstream; prion protein 2 (dublet)		

For research use only.

GENE DESCRIPTION	prion protein 2 (dublet)
GENE ID	23627
ORIGIN SPECIES	Human
ACCESSIONS	BC043644
VECTOR	pDONR223
NCBI SYMBOL	PRND
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggagcaccctgagctgggtggctggccactgtctgcatgctgtcttcagccacctctctgcggtccagacgaggggcatcaagcacagaatcaagtgaaccggaaggccctgccagcactgccagatcactgagccagggtggtgagaaccgcccgggagccttc
atcaagcaaggccgaagctcgacattgacttggagccgagggcaacaggtactacgagggcaactactggcagttcccgatggcaccactacaacggctgctctgaggctaatgtgaccaaggaggcattgtcaccggctgcatcaatgccaccaggcggcgaaccagggg
gagttccagaagccagacaacaagctccaccagcagggtgctctggcggctgtccaggagctctgtccctcaagcattgcgagttttggttgagagggggcgaggacttcgggtcaccatgcaccagccagtgctcctctgccttggcttggatctggtcatggtgaataa