



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

IL18BP-interleukin 18 binding protein Gene

Reference: PTXBC044215

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

IL18BPA; interleukin-18-binding protein; IL-18BP; MC51L-53L-54L homolog gene product; tadekinig-alfa; interleukin 18 binding protein

For research use only.

GENE DESCRIPTION	interleukin 18 binding protein
GENE ID	10068
ORIGIN SPECIES	Human
ACCESSIONS	BC044215
VECTOR	pDONR223
NCBI SYMBOL	IL18BP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccatgagacacaactggacaccagacctcagcccttctgggtctctctctctgtgtgccacgtctgctactctctctgtgtagacccacactgtctctcgagaccaccacagctgccactccctcagttagaagcacaaaggaccctgccctcccaagccccagtgftccacag
ctaagcagtggtccagcattggaagtgacctggccagaggttggaagtggccactgaatggaaacgtgagcttatctgtgtggcctcgacggccttcccacactcagcattccttactggtctgggcaattgttctcattgagcacctccaggccgactgtgggagggggagcaccagccg
ggaaacgtggggagacagctatgcagctctgtgcaagccctgtgtctggagcagctgacccctcgacacactcaacacttctctgtctgtgacccctgaacaggtgttcacagctgcacgtctgccctgcccaactctgtgggtgggtgagggaacacttgcctccaccacaag
aacctctccctccctcacacacagcagctctcagacacacagggttaa