

ANKRD13D-ankyrin repeat domain 13 family, member D Gene

Reference: PTXBC044239

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ankyrin repeat domain-containing protein 13D; ankyrin repeat domain 13 family, member D; ankyrin repeat domain 13D		

For research use only.

GENE DESCRIPTION	ankyrin repeat domain 13 family, member D
GENE ID	338692
ORIGIN SPECIES	Human
ACCESSIONS	BC044239
VECTOR	pDONR223
NCBI SYMBOL	ANKRD13D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcctgtggtggcgtggcaggttcaaggcaactgtggtgagtgagagacccgctctccctgggtgaccaggtgacccccatcatcgacctaagtcgcatcagcaacgtcactttgccaagctgcgcgactcatcactctgcgccttcacactggcttcccggtcaaaattga
gattcccttttccacgtgctcaatgcccgcacaccttcagcaacctgtgtggtgtgagagccctgagctccgtgtggtggtgcggccccccagctctgctgctgcgcgcatcagggaacctttcccggtgcgaggtggaccccacccgtgttgaagtgcccaacgggtacagcgtgctgg
gcatggagcgcaacgagccccctcgggacgaggacgatgacctcctgcagttcgccatccagcagagcctgcttgaagcgggactgagggcgagcaggtgaccgtctgggaagccctgaccaacacccggcccggtgccgcccctctcccagccacgggttatgaggaaca
gcttcagctggagggccctccaggaaagcctgcagctgtccacagagccaggggcccgagatccctccagacaccccgccccgggtccacccagcttgaagagcagctgcgcctggccctggagttgtcttcaggagcaggagagcgggagcggcgccggca
gcaggaggaggagacttacagcgatcctgcagctgtcactcactgagcactga