

## MUC20-mucin 20, cell surface associated Gene

Reference: PTXBC044243

|                               |                                                                                 |                     |                |
|-------------------------------|---------------------------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | MUC-20; mucin-20; transmembrane mucin MUC20S; mucin 20, cell surface associated |                     |                |
| <b>For research use only.</b> |                                                                                 |                     |                |

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | mucin 20, cell surface associated                                                                |
| <b>GENE ID</b>          | 200958                                                                                           |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC044243                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | MUC20                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgggctgtctctggggtctggcctctgccctttctctctgtctgggaggttggggtctctgggagctctgcaggccccagcaccgagagcagacactgcgatgacaacggacgacacagaagtgtccgctatgactctagcaccgggcccagccgctctggaactcaaacgctg  
agcgtgagacctctcttagggcctcaacccagccggccccattccagaagcagagaccaggggagccaagaatgtccctgcaagagagaccaggagtttcacaaaaacatctcccaactcatgggtctgctgacccacctcgtggagacatcagccgccagtgccagccccg  
aggagctgggaatgaccacagttcagaccatcacaggcagtgatccggaggaagccatcttgcacacctttgcaccgatgacagctctgaagaggcaaaagacactcacaatggacatattgacattggctcacacctccacagaagctaaaggcctgtcctcagagagcagtgctcttc  
cgacggcccccatccagtcacccccgtcacgggctcagagagcagcgccctctccgacggcccccatccagtcacccccgtcacccccgtcatgggtcccccggatctgatgtcactctctcgtctg  
aagccctgggtgactgtcacaacatcgaggttaattgagcgtacacagaataagaaacaaactccagcatccctggggcctcagacatagatctcatcccaaggaggggtgaaggcctcgtccacctccgacccagctctgctgactccactgaagaaaaaacacat  
cactgaggtcacagcctctgcccagacctgtccacagccggcaccacagagtcagctgcaacctatgccacggttgggacccccactccccactaacagcgcacacagaagaagtgacagcaccggggccagcaccctcagtgagctctgtgtcagtagcaggaatccct  
ggaagaaactcagccctctctgttgagacaccaagttacgtcaaaagtcagagcagctccggtctccatagaggctgggtcagcagtgaggcaaaacactctcttctgtgggagctctgcttctctacagccccctgggaagccgcccccaagaacttcacccctcagagacaccg  
accatggacatcgcaaccaaggggccccctccaccagcaggacccctctctctgtctctccgactacaaccaacagcagccgaggagacagacacacttagccaagatcacaacctcagcgaagaccacgatgaagcccccaacagccagccccagctgcccggacgag  
gccgaccacagacgtgagtgacaggtgaaatgaggtttctctctctgctgctgagtggtgtcccccgaagacctcactgacccagagtggcagaaaggctgatgcagcagctccaccgggaactccagccccagcgcctcacttcaggtctcttactgctgtcaggagag  
gctaa