

TSPAN33-tetraspanin 33 Gene

Reference: PTXBC044244

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PEN; PEN; TSPAN-33; tetraspanin-33; penumbra; proerythroblast new membrane; tetraspanin 33		

For research use only.

GENE DESCRIPTION	tetraspanin 33
GENE ID	340348
ORIGIN SPECIES	Human
ACCESSIONS	BC044244
VECTOR	pDONR223
NCBI SYMBOL	TSPAN33
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaggagacccggcgccgcccgcgggaggttctcttcagccccgctggtgaatacctgctcttcttcaacatgctcttctgggtgatttccatggtgatgggtggtggtgtctacgctcggctaataagcatgcagaagcagccctagcctggcga
gtggaccctgccatctgctgacgtgtgtgggtgtcctcatgttctgctcacccttctgtgctgcatgggtccctccgcgagaacatctgctcctcagacgttctcctctgctcaccgctgtgttctgctgagctggcgctgggacccctggcttctctcagacaaggctcga
gggaaagtgagtgagatcaacaatgccattgtgcactaccgagatgacttggatctgcagaacctcattgatttggccagaaaaagtttagctgctgtggagggtattctacaaggactggtctcagaacatgtatttcaactgctcagaagacaacccagtcgagagcgtgctctgt
gccttactctgttctgcttgcctactcctgaccaggcagtgatcaacactatgtgtggccaaggtatgcaggcctttgactacttggaaagctagcaaaagtcacacccaatggctgtattgacaagttggtaactggatacagcagaacatttcttacttgggtgtggtcttaggcctggcc
atccccagctgtgtgggaattctgttcccagatcctagtgaatcagatcaaaagatcagatcaagctacagctctacaaccagcagcaccgggctgacctggtactga