



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

GPM6A-glycoprotein M6A Gene

Reference: PTXBC044612

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

GPM6; M6A; neuronal membrane glycoprotein M6-a; glycoprotein M6A

For research use only.

GENE DESCRIPTION	glycoprotein M6A
GENE ID	2823
ORIGIN SPECIES	Human
ACCESSIONS	BC044612
VECTOR	pDONR223
NCBI SYMBOL	GPM6A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagagaataatggaagaggacacacaaaaagggtgttttgaatgctgtatcaaatgcctggggggcattccctatgctctctgattgccaccatctgctctatgcgggtgtgcccctgttctggctcgggctatgaagcgcttctggaaatgtcaacattctgcaaacctattt
gagatggcaagaactcgtgagacacactggatgtttttaccatgattgacatctttaatagtatgctgaltacggatcgcagctgcgtcttftgtatggcatttgcgtatgtggaaggtttcttcacaactgggggcacaaagatctctatgggaattcaaaatcacacttctggcgcatgt
gtgagcgtttgttcattatgctgacatattttcatgttgcctggctggagctcagcgtttcaactcactgcgattttacatgtactcaatctgtggacatctgcgggaacaccacattatggagggagcaaaactctctgtgggaacttctgatttggaaattgacaatttggagga
aaagaaaatttactgtctgacgaalltctgagatgtgcgaattcactgactgagctgaacatcagcttccacttctttatgtggcattctggagatggcgagtcagatgctatgttgcactactatattgttctgtctccaactggccctatgtgaagacgcctgcgggatcgagaat
atgaagacatcaagtctcgaaggaagaacgaagatgcagatcactctatcctcctcaaaaggcgatcactacacataa

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1