

C10orf93-chromosome 10 open reading frame 93 Gene

Reference: PTXBC044661

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TPR repeat-containing protein C10orf93; C10orf93; C10orf123; C10orf124; C10orf92; TTC40; bA288G11.4; bA288G11.5; bB137A17.2; bB137A17.3; cilia- and flagella-associated protein 46; protein CFAP46; tetratricopeptide repeat domain 40; tetratricopeptide repea		

For research use only.

GENE DESCRIPTION	chromosome 10 open reading frame 93
GENE ID	255352
ORIGIN SPECIES	Human
ACCESSIONS	BC044661
VECTOR	pDONR223
NCBI SYMBOL	C10orf93
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacctggtcatcagcaggagctggcccgcgcgagagccagcaagatctgctccttgaagaaggcctacaggttgatcaaatcgccaacctagggaatcggaagttgacccctcagagagctcagccagacctgtttgtctgtgtgcagagcagccctgaagatgag
gcagccagagggtgagcaggactgcatccaaatgtactcaagggtgaaggcgcacatcaccagtttctggccgagcgcacctgtgcagggccagatgtgtgccccgaagtcgagcaaaaacctggagggaatttgaaattgctgactgagtacatgaaggccataaactttgcca
aaggagaaccgaggtactactttttgtgtacaatgcatcagtccttactgagcagatgggtgagccgtttcctcaagcctggaatcgtcaccatctgacccagccttcccaaatcataaacgtgctgagtcagactgaggagaagacaaggagtgccgtgctgagctgatgctggaac
ttctggagtgttatctgcaagccggaagaaaggagaggtccaggttctgtccacggcagctcgttcaatgaatctcagctgccacagaataccggcagatattctctgtatgttaaatcagcgtattagaataagagacgtaagaagtaattttctaaagtcattgttcagaaaact
gcgtagaaattacataa