

C10orf96-chromosome 10 open reading frame 96 Gene

Reference: PTXBC044830

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0628 protein C10orf96; C10orf96; coiled-coil domain-containing protein 172; coiled-coil domain containing 172		

For research use only.

GENE DESCRIPTION	chromosome 10 open reading frame 96
GENE ID	374355
ORIGIN SPECIES	Human
ACCESSIONS	BC044830
VECTOR	pDONR223
NCBI SYMBOL	C10orf96
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcttggagtcctgtttcagcacatcatctcaccgagcatcaggcggaggagagtcgctgttgatgcgagaagtaaggctggaaataaccagatgctgtgaaaaaataagaaagcaacggaggagctgaatgaagagaaatcaagctggaatctaaggtcaacagttttga
aaaatcctctcttacagcttttgaagctcatgaaatgctttgaaaaacagtagagtgaaattacaacataggaatgcttctcaaacctttgaggtatataagaaacaaatgtagaggaggaagacaatttataaggaaattacagactttaaatgattgaaatacaagaa
aagagagcttttgatgaagaaatgtcaagattgaaatctgacttagaaaaccaagcaaacatgttgaagaatgaagtaagtcgaatgatagtagccagtttaataaactcaaaaacaaaagagtgaaatgatacaagaattatttactctccaaagaaaactaaagtgttgaa
gatgaagagaatgaatccattgtactacaaatctagaggcagaaaaataaaatcagtgaaaagcctcaaatgatacagaatgcttaagacttaaaaaagaattagaactttataagggaagatgacatggaagtgtttatgaagctctccaaacagaaatagaattttggaggtgaca
ttggcacagaaagatcttcaggaagcaataa