



Reference: PTXBC045113

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	leucine-rich repeat transmembrane neuronal protein 1; leucine-rich repeat transmembrane neuronal 1 protein; leucine rich repeat transmembrane neuronal 1		

For research use only.

GENE DESCRIPTION	leucine rich repeat transmembrane neuronal 1
GENE ID	347730
ORIGIN SPECIES	Human
ACCESSIONS	BC045113
VECTOR	pDONR223
NCBI SYMBOL	LRRTM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]