

SLAIN1-SLAIN motif family, member 1 Gene

Reference: PTXBC045177

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C13orf32; SLAIN motif-containing protein 1; SLAIN motif family member 1		
For research use only.			

GENE DESCRIPTION	SLAIN motif family, member 1
GENE ID	122060
ORIGIN SPECIES	Human
ACCESSIONS	BC045177
VECTOR	pDONR223
NCBI SYMBOL	SLAIN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggataataattacaggacctcactgatgttcagatcatggctcgtcgaagaagaaagtcaggcaagattatgcttctactcagcatctgtatcaagacatagttccagtggtcattgagttcaggaaaaaaggacatgtagtatgaagaatatgaccaatacagctcggaggat
gaagaggaatttgatcatttgccaccacctcagcctcgttccaagatgttcccttccaagagggaattcccatcagacagcttctccagcattcgggagtgtagggaggagccccagttccagatatttccctcaataaattaccagcagcaacagatatttcacctcaagccaaactcc
agatcagcaaccaaataggaccaatggagataagctccgaagaagtatgcttaacctagcccgatgccaagtacaactgccattagtagcaacattagttctccgggtcaccgtgcgaataatgcagagttttgactcaagcttgcattggagctggaaatggaattcaagaatacattgt
attccatcacgggacagcttcaacacaggggtccacagcgtggggcatttccagtgctatccgacagcctttaaagccacagcctatgtgagtcacacgttcaaggcagcagtaacatgcctttatcaaacggccttacagctgtattccaacaggaatccccacccaacgaacaaagct
gcagcttctgggataatgggtcgcagtcgactccaagacctcgttggcaataaatgggagtaacctgcctcgaagcaaaatgcacaacctgttagaagttttctcagcctccaaagcctctgtcttactcagcactctgaggatggaaattggagagatggttctactaa