

AFG3L1-AFG3 ATPase family gene 3-like 1 (*S. cerevisiae*) Gene

Reference: PTXBC045753

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	1700047G05Rik; 3110061K15Rik; AFG3-like protein 1; AFG3(ATPase family gene 3)-like 1; AFG3-like AAA ATPase 1		

For research use only.

GENE DESCRIPTION	AFG3 ATPase family gene 3-like 1 (<i>S. cerevisiae</i>)
GENE ID	172
ORIGIN SPECIES	Human
ACCESSIONS	BC045753
VECTOR	pDONR223
NCBI SYMBOL	AFG3L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgagtcctggaggagacggaggcaggagaggagggaagcaggataacgttgccctggtggaggcggatgcagaaggggggactccctgggatgacgaggatttccgcagctctggccctttgggggcaggcgttgccatgggattttctacctctatttcgagatcctggaaga
gaaatcacgtgggaagcattgtacagtattacctggccagaggtctgtggaccggctggaagtcgtgaacaacaatctgtgcgtgttattcctgcccctgggacctctctgaagtgagggtgaatttaaggctgagtattgcagacataagtttatttcgtgtaaaatgtagttttctatttt
tttcagtaa