

ZNF283-zinc finger protein 283 Gene

Reference: PTXBC045755

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HZF19; HZF41; zinc finger protein 283; zinc finger protein 41; zinc finger protein HZF19		

For research use only.

GENE DESCRIPTION	zinc finger protein 283
GENE ID	284349
ORIGIN SPECIES	Human
ACCESSIONS	BC045755
VECTOR	pDONR223
NCBI SYMBOL	ZNF283
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcgagagctgttctggtttcttggaattctgtgcttcaccaatagaggaatccatggagcattaataattcttgcattccagaacctgactgatgggttggtgacattcagggatgtggccatcgacttctctcaggaggagtgggaatgcctggaccctgctcagagggtgtacgtggatgtaattgtggagaactatagtaacttggtgtcactgggtaaggatcatctgcctgaaataa