

LOC374491-TPTE and PTEN homologous inositol lipid phosphatase pseudogene Gene

Reference: PTXBC045761

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	TPTE and PTEN homologous inositol lipid phosphatase pseudogene
GENE ID	374491
ORIGIN SPECIES	Human
ACCESSIONS	BC045761
VECTOR	pDONR223
NCBI SYMBOL	LOC374491
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaaacgtgctactgagatatgttgatatttgcacaagtgaacatctctacaactggaatctccccaagacggatactctttataaagattcattatttattcgattcggtgttggaacagcgcatgtatgtgatctaaaattccaaatagtaattggagaaaaagatactgcatacaca
ttgaaacagctgggtgtattaactaatgtatatgacagtcctctctgtacgatgtgaaagtgcagttttctcttcaaatcttccataatactatgacaattgttcattttctctctggttcaacacatctttattcaaaataacagagagaatttagatatatttataaacttgaagtcattaa