

LVRN-laeverin Gene

Reference: PTXBC045809

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APQ; AQPEP; TAQPEP; aminopeptidase Q; AP-Q; CHL2 antigen; laeverin		

For research use only.

GENE DESCRIPTION	laeverin
GENE ID	206338
ORIGIN SPECIES	Human
ACCESSIONS	BC045809
VECTOR	pDONR223
NCBI SYMBOL	LVRN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggtggaaatttcaaaacaagtgaatacagggaactctttgacatatttacttacagcaaggagcgtctatggcccgatgcttctgtttctgaatgagcatttattgtcagtcgactcaagtcataattgaagacattttctactcaaacgctgagcaagatgatctatggagcattt
tcaaatggccatagatgaccagatgacagttattttgccagcaacaataaaaaacataatggacagttggacacaccagagtggtttccagtgatcactttaaatgtgtctactggcgtcatgaaacaggagccattttatcttgaacattaaaaatcggactcttaccagcaatgacaca
tggattgtccctattctttggataaaaaatggaactacacaacctttagtctggctagatcaaagcagcaaatattccgaaatgcaagttcagattctgacatgactgggtgattttgaattgaaatgactggatattatagagttaattatgataaattagggttggagaaactaaatcaaca
acttgaagagatcctaagatgagataa