

OCC-1-overexpressed in colon carcinoma-1 Gene

Reference: PTXBC045812

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OCC-1; overexpressed in colon carcinoma 1 protein homolog; overexpressed in colon carcinoma 1 homolog; chromosome 12 open reading frame 75		

For research use only.

GENE DESCRIPTION	overexpressed in colon carcinoma-1
GENE ID	387882
ORIGIN SPECIES	Human
ACCESSIONS	BC045812
VECTOR	pDONR223
NCBI SYMBOL	OCC-1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactcaagaatcaagagcttgctcatcagtttgaagggaattggctccgtgggacgttgtaatgtgcacagacattccaaggaaattctaacagtcacccttcctttgctatcccccaatcttaagtgatataaaaacccctgggtacatatgttggttaataagaaggaattggtaaacagtacacttggttatggaactttctgtggccacctacgaaagacaagttaa