

VBP1-von Hippel-Lindau binding protein 1 Gene

Reference: PTXBC046094

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HIBBJ46; PFD3; PFDN3; VBP-1; prefoldin subunit 3; von Hippel-Lindau binding protein 1; VHL binding protein 1		

For research use only.

GENE DESCRIPTION	von Hippel-Lindau binding protein 1
GENE ID	7411
ORIGIN SPECIES	Human
ACCESSIONS	BC046094
VECTOR	pDONR223
NCBI SYMBOL	VBP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccgttaaggacagttgtggcaaggagaaatggccacaggggaatggcgccgctccactggggattcctgagggcgtgtttgtggaagatgtagattcctcatgaacagcctgggaatgagactgcagatacagtattaaagaagctggatgaacagtagcagaagtat
aagtttatggaaactcaacctgtctcaaaagaaagaaggtctaaaggtcagattcctgaaattaaacagactttgaaattctaaatacatgcagaagaaaaagagtcaccaactcaatggagaccagattctgtgagataacctgtattgcaaagcttcagttctctaccgataa
agtggtctgtgttggggcctaattgaatgctgaatgatattgatgaagctcaggcattgttgaaaagaattatcgactgcccaagaatcttgattcctggaggaagacctgactttctcgagatcaatttaccacagaagtcaatatgccagggtttataattgggatgtaaa
aagaagaacaaggatgactctaccaagaacaaagcataa