



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

WFDC2-WAP four-disulfide core domain 2 Gene

Reference: PTXBC046106

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EDDM4; HE4; WAP5; dJ461P17.6; WAP four-disulfide core domain protein 2; WAP domain containing protein HE4-V4; epididymal protein 4; epididymal secretory protein E4; epididymis-specific, whey-acidic protein type, four-disulfide core; major epididymis-speci		

For research use only.

GENE DESCRIPTION	WAP four-disulfide core domain 2
GENE ID	10406
ORIGIN SPECIES	Human
ACCESSIONS	BC046106
VECTOR	pDONR223
NCBI SYMBOL	WFDC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctctgtgttgctgcttagggccgctagccgcccctctccctcagctctgtctgttgcgtttaccctagtctcaggcacaggagcagagaagactggcgtgtgtcccgagctccaggctgaccacaactgcacgcgaagtagtgcgtctcggacagcgaatgccgcgacaacctcaa
gtctcgtcgtcgtcgtgctgtccactcttctgtctctgcgtcccaatgataaggagggttctgtccccaggctgtgaacatlaattcttccccagctcgcgcctctgtcggaccacgtgccagggtggacagccagtgtcttccggccagatgaaatgtctccgcaatggctgtggggaagggtctctgtgtca
ctcccaatttctga