

TMEM184C-transmembrane protein 184C Gene

Reference: PTXBC046128

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TMEM34; transmembrane protein 184C; transmembrane protein 34		

For research use only.

GENE DESCRIPTION	transmembrane protein 184C
GENE ID	55751
ORIGIN SPECIES	Human
ACCESSIONS	BC046128
VECTOR	pDONR223
NCBI SYMBOL	TMEM184C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctgcactgtacctggaggaaactggagacagtggaatcgacctttagtagcgggtcatctacgtggtgcaatagtggttgcgggtccctatgcgtgtgggaattacagaaactggaggttgaatacacaccaaggcttggttattgctggaatcttttgcgttgactattcctatatca
ctgtgggtgataattgcaacacttagtcattatacacacactgaactacaaaaccaataaaggattcttggatggtacatttacagtttagatagttggatagcttgaatatcccggaattgcaatatatgtgatactgcagagaatgctatgaagcttatgtaattacaactttatggg
attccttaccattatctaaactaaccgggtatccaaatctgggtattaatccttgaagccaagatcaacagaacatttccctctttagttgctgtcccatgggctatgggagagattgctgttttaggtgcaactagggtgtattacagtacacagttgtcagaccttcaccaccatcgttgcttt
aatctgtgagctcgttggtatatatgacgaagggaactttagcttttcaaatgcttgacttattggttataataacaacatgtcacagttgttgcctgtattgtctcgtctttataaaagtactaaaagaagaactgagccaatccaacctgttgcaaatctttgtgtaagctggtggttt
tgtttctttttggccttaa