

COMMD2-COMM domain containing 2 Gene

Reference: PTXBC046131

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC042; COMM domain-containing protein 2; COMM domain containing 2		
For research use only.			

GENE DESCRIPTION	COMM domain containing 2
GENE ID	51122
ORIGIN SPECIES	Human
ACCESSIONS	BC046131
VECTOR	pDONR223
NCBI SYMBOL	COMMD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgctggaattgtccgaggagcataaggaacacctggccttcctgctcaagtgagacgcgggtgtgctgctggaattctgagacgcggcgcaaaccccaaaatctacgaaggcgccgagaaaactcaatgtgagtagtgacactgtccagcat
gggtgtggaaggatlaacgtatctcctcactgagagctcaagctcatgattctgaactggatttccaagactctgtttgtgtctgggattctctgaagaattaaacaaatgtgtctcagctttatctggacaacagaaaagagatcagaacgtctctgagtgaaatggcaccagcctccagtt
atcataacctgaatggcgactagatgtacagcttgcaagtagaagctcagggaacagattaaaccagcagtgactataaagctacaccttaatacaaaatggagatcacacaccaaagtctctgagacagaccagccacctgtccatttgggttcaacaactggaacaagcattggaag
agatgaagacaaaccactgtaggagagttgttcgcaacatcaagtag