

NDUFS6-NADH dehydrogenase (ubiquinone) Fe-S protein 6, 13kDa (NADH-coenzyme Q reductase) Gene

Reference: PTXBC046155

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NADH:ubiquinone oxidoreductase NDUFS6 subunit; CI-13kA; CI-13kD-A; CI13KDA; NADH dehydrogenase [ubiquinone] iron-sulfur protein 6, mitochondrial; NADH dehydrogenase (ubiquinone) Fe-S protein 6, 13kDa (NADH-coenzyme Q reductase); NADH-ubiquinone oxidoreduc		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) Fe-S protein 6, 13kDa (NADH-coenzyme Q reductase)
GENE ID	4726
ORIGIN SPECIES	Human
ACCESSIONS	BC046155
VECTOR	pDONR223
NCBI SYMBOL	NDUFS6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgccgatgacctctgccgctgctgaaccgtgtggcgagggcgccggagcctgccctggcgccaggtgttcgggggtgcgggtctcgcgcgaccggggagaggtcacgcacactggccaggttatgatgataaagactacaggagaattcgtttgtagtcgtc
agaaagaggtgaatgaaaacttgcattgattgatagcagagcagccctgagcgaagggtggagactcgggtgatagcgtgcgatggcgcgccgggagctcttggccacccaaaagtgtatataaacttgacaaagaacaaaaccggcacatgcggttactgtggctccagttc
agacagcaccaccactag