

DIP2A-DIP2 disco-interacting protein 2 homolog A (Drosophila) Gene

Reference: PTXBC046176

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf106; DIP2; disco-interacting protein 2 homolog A; DIP2 disco-interacting protein 2 homolog A; DIP2 homolog A; disco-interacting protein 2A; disco interacting protein 2 homolog A		

For research use only.

GENE DESCRIPTION	DIP2 disco-interacting protein 2 homolog A (Drosophila)
GENE ID	23181
ORIGIN SPECIES	Human
ACCESSIONS	BC046176
VECTOR	pDONR223
NCBI SYMBOL	DIP2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgggtaccacccatcgacattgagacctgtgtcatccgagcacacaggagcatcgctgagtggtgccgtattcacctggaccaacctgctgggtggtggtggtgagctggatgggctagagcaggatgccctggacctgggtggcctggtgaccaacgtggtgctggaggagcacta
cctggctgtgggagtggtggtcatcggtgaccaggggtgatccctatcaactctcgggtgagaagcagcgcacgtcacctgcgggacggccttcctgctgaccagctggacccatctatgtcgctacaacatgtga