

OVOS2-ovostatin 2 Gene

Reference: PTXBC046190

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ovostatin homolog 2; ovostatin 2		
For research use only.			

GENE DESCRIPTION	ovostatin 2
GENE ID	144203
ORIGIN SPECIES	Human
ACCESSIONS	BC046190
VECTOR	pDONR223
NCBI SYMBOL	OVOS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggttatagatgtaaaatgctatcaggatttactccaacatgtcatccattgaagagcttgaacaaggcccaagtgatgaagactgaagtcaagaatgaccatgttcttttctacttggaaaatgttttggcgagcagacagtttcttttctgttgagcagagcaacctgtgttca
acattcagccagccccagcatggtctacgattactacgaaaaagaagaatagccctagctttttaccacatcaacagtagttcagttccgagtga