

IFIH1-interferon induced with helicase C domain 1 Gene

Reference: PTXBC046208

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AGS7; Hlcl; IDDM19; MDA-5; MDA5; RLR-2; SGMRT1; interferon-induced helicase C domain-containing protein 1; CADM-140 autoantigen; DEAD/H (Asp-Glu-Ala-Asp/His) box polypeptide; RIG-I-like receptor 2; RNA helicase-DEAD box protein 116; clinically amyopathic		

For research use only.

GENE DESCRIPTION	interferon induced with helicase C domain 1
GENE ID	64135
ORIGIN SPECIES	Human
ACCESSIONS	BC046208
VECTOR	pDONR223
NCBI SYMBOL	IFIH1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgaatgggtattccacagacgagaatttcgctatctcctgcttcaggccagggtgaaaatgtacatccaggtggagcctgtgctggactacctgaccttctgcctgcagaggtgaaggagcagattcagaggacagtcgccacctccgggaacatgcaggcagttgaact
gctgctgagcaccttgagaaggagctcggcaccttggtggactcgggaattcgtggagccctccggagaaccggcagccctctggccgccgctacatgaacctgagctcacggacttgcctctccatcgtttgagaacgctcatgatgaatctccaactgctgaacctcttc
agccactctggtggacaagcttctagtagacgtcttgataagtgcattgagggaggaactgttgacaattgaagacagaaacggattgctgctgcagaaaaaatggaaatgaatcaggtgtgaagagagctactaaaaaggattgtgcagaaagaaactggttctctgcatcttga
atgttctctcaaacaggaacaatgaactgtccaagagtaacaggctctgattgctcagaagcaatgcaggtatttgaatttactgagggaagattctcaattctgcctag