

RPAIN-RPA interacting protein Gene

Reference: PTXBC046349

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HRIP; RIP; RPA-interacting protein; RAP interaction protein; nuclear transporter; RPA interacting protein		

For research use only.

GENE DESCRIPTION	RPA interacting protein
GENE ID	84268
ORIGIN SPECIES	Human
ACCESSIONS	BC046349
VECTOR	pDONR223
NCBI SYMBOL	RPAIN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggagtcgttgaggtcctccgcgcctcctgtacaaactgtggctcgccttgaaagagccttccgagagatgcctggagagaatgagaacagccggacagctcctaaacaggtaccgagctggaagcagtgaggccagggaattctcagaacagcttctagttcaagagtgatggaagaagagtggaatccttgcagtcagtgagaaattgccagaagactggctcagctgaggagctgatagacatggctgtgctggagaaattcaacaggagctgatcaaccaaggcctgtga