

TMEM37-transmembrane protein 37 Gene

Reference: PTXBC046362

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PR1; voltage-dependent calcium channel gamma-like subunit; neuronal voltage-gated calcium channel gamma-like subunit; voltage-dependent calcium channel gamma subunit-like protein; transmembrane protein 37		

For research use only.

GENE DESCRIPTION	transmembrane protein 37
GENE ID	140738
ORIGIN SPECIES	Human
ACCESSIONS	BC046362
VECTOR	pDONR223
NCBI SYMBOL	TMEM37
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctgtacgcagcgtggcgcccttgccgtgtggccgccaattttggcctggagttcctcatggtgccagttgtgcgagacaacactcacagtgcagtggtgatgggtccatcctcctctgtgtctttcctcctcctcgggggctcctgggtttgatcctcctc
aggaaaccaagtacacatcatcggcttaccctaattgtttgtgcgaattcactgcctcctcctcctctgaacgccatcagcggccttcaacaacagcagcaccatccctgggaatga