

CPNE8-copine VIII Gene

Reference: PTXBC046366

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	copine-8; copine VIII; copine 8		

For research use only.

GENE DESCRIPTION	copine VIII
GENE ID	144402
ORIGIN SPECIES	Human
ACCESSIONS	BC046366
VECTOR	pDONR223
NCBI SYMBOL	CPNE8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatccttaccactgaatgcctatggtatggcactaaaagcagtgaggagaaattgtcaagattatgacagtataaaatgttccagctctaggatttgggtgcaaaactgcctccagatggaaggatatctcagaatttgccttgaatgggaatcctcaaaaccctactgtgatggcattg
agggggtcatggaggcttattacaggagctgaaatctgtacaactatatgggccccaactttgtcctgttaataatcatgtacgaagatatgctctctgttaaaggatggctccagatatttgtgcttctgattgttacagatgggtgttatctcagatatggccagactaaggagtcctagat
taatgcctcaaaacttccaatgtcaataattatagtaggtgtggaccagcagaatttgcgaatggtcgaattggaatggagatgatgaagagctcctctagaggaaaatagctgaaagagacattgtgcagtttggccattcagggaattatattgacagaagtggaaaccatactgagc
atggctagattggcctaaagatgtcctagctgagatccctgagcagtttctctctatagagcccagggaatcaagccatcacctgcgcctccccatacaccacacacatgtgttacagactcaaatatga