

## TMEM77-transmembrane protein 77 Gene

Reference: PTXBC047025

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TMEM77; CORD21; PRO180; WWFQ154; DNA damage-regulated autophagy modulator protein 2; RP5-1180E21.1; damage regulated autophagy modulator 2; transmembrane protein 77; DNA damage regulated autophagy modulator 2		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	transmembrane protein 77
<b>GENE ID</b>	128338
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC047025
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	TMEM77
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

atgctaaatattcggcagtttatgcattctaccatttatgttcgttataagcaagttcatgctctgagctctgaagagaagcttatcatcaaataaacaaggctggcctgtactggaatactgagttgttaggacttctattgtggcaactccagaaacaacccttttctgcacatgtaa  
gtggagctgtgcttaccttggatfaggctcattatataatgttttcagaccatccctccacaaaatcagcccaaaatccatggcaacaagctctctggatcagactgtgttggttatctggtggagtaagtcacttagcatgctgactgctcatcagtttgcacagtggaatttggg  
actgattagaacagaaactccatggaaacccgaggacaagggttatgtgcttcacatgatcactactgcagcagaatggtctatgtcatttctcttgggtttttctgacttacctctgatttcagaaaaatctttacgggtggaagccaattacatggattaaccctctatgacactgcac  
cttgccttatacaatgaacgaacacggctacttccagagatattga