

IL33-interleukin 33 Gene

Reference: PTXBC047085

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf26; DVS27; IL1F11; NF-HEV; NFEHEV; interleukin-33; DVS27-related protein; interleukin-1 family member 11; nuclear factor for high endothelial venules; nuclear factor from high endothelial venules; interleukin 33		

For research use only.

GENE DESCRIPTION	interleukin 33
GENE ID	90865
ORIGIN SPECIES	Human
ACCESSIONS	BC047085
VECTOR	pDONR223
NCBI SYMBOL	IL33
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagcctaaatgaagtattcaaccaacaaattccacagcaaaagtgaagaacacagcaagcaaaagcctgtgttcaagctgggaaatccaacagaaggccaaagaagttgccccatgtactttatgaagctccgctctggccttatgataaaaaaggagcctgttacttagg
agagaaaccaccaaagccttcactgaaaaaggtagaagcacaagaacatctgtgtactcgtcgtgtcaacagcagtcactgtgagtgcttggcctttgtatcaggggtccagaaatatactagacacttcattcaagatcacaggaaattcactattacagatctt
gcttcttaagcacatacaatgatcaatccattcttggcttggaggtgaaagttagagatatgtgaagactgaaaaaagatgaaagaagataaggtgttactgagttactatgagctcaacacccctcaatgaatcagggtgacggtgttggatgtaaggttaagtgaaccctg
agtcctacaaaagacttctgttgcacccaacaaggaacactctgtggagctccataaggtgtgaaaaaccactgccagaccaggcctcttctcctcataatatgcactccaactgtgttcatgtgaatgcaagactgacctggaggtttataggtgtaaggataatcattctctg
attaaagtagactctctgagaatttgtgtactgaaaatatctgtttaagctctctgaaacttag