



Reference: PTXBC047129

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	dimethylaniline monooxygenase [N-oxide-forming] 1; FMO 1; Flavin-containing monooxygenase 1 (fetal liver); dimethylaniline oxidase 1; fetal hepatic flavin-containing monooxygenase 1; flavin containing monooxygenase 1		

For research use only.

GENE DESCRIPTION	flavin containing monooxygenase 1
GENE ID	2326
ORIGIN SPECIES	Human
ACCESSIONS	BC047129
VECTOR	pDONR223
NCBI SYMBOL	FMO1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr