

EMP1-epithelial membrane protein 1 Gene

Reference: PTXBC047300

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CL-20; EMP-1; TMP; epithelial membrane protein 1; tumor-associated membrane protein		

For research use only.

GENE DESCRIPTION	epithelial membrane protein 1
GENE ID	2012
ORIGIN SPECIES	Human
ACCESSIONS	BC047300
VECTOR	pDONR223
NCBI SYMBOL	EMP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttggtattgctggctgtatctttgtggccacatcgctactgtttattatgctattttagcaccattgccaatgtctgtgtgttccaatacggtagatgcatcagtaggtcttttgaaaaactgtaccaacattagctgcagtgacagcctgtcatatgccagtgaagatgccctcaagacagt
gcaggccttcattgattctctatcatctctgtgtcattgccctcctgtgtcttcctgtgttccagctcttcacatggagaagggaaccgggttctctctcaggggccaccacactggtgtgtgctgtgtgcattcttgggggtgtccatctacactagtcattatgcgaatcgtgatggaacgca
gtatcaccacggctattcctacatcctgggctggatctgtctctcagcttcacatcggtctctctatctgttctctgagaaagaaataa