



Reference: PTXBC047318

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSS; MPS3A; SFMD; N-sulphoglucosamine sulphohydrolase; heparan sulfate sulfatase; mucopolysaccharidosis type IIIA; sulfoglucosamine sulfamidase; sulphamidase; N-sulfoglucosamine sulfohydrolase		

For research use only.

GENE DESCRIPTION	N-sulfoglucosamine sulfohydrolase
GENE ID	6448
ORIGIN SPECIES	Human
ACCESSIONS	BC047318
VECTOR	pDONR223
NCBI SYMBOL	SGSH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr