

RCHY1-ring finger and CHY zinc finger domain containing 1 Gene

Reference: PTXBC047393

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARNIP; CHIMP; PIRH2; PRO1996; RNF199; ZCHY; ZNF363; RING finger and CHY zinc finger domain-containing protein 1; CH-rich interacting match with PLAG1; E3 ubiquitin-protein ligase Pirh2; RING finger protein 199; androgen-receptor N-terminal-interacting pro		

For research use only.

GENE DESCRIPTION	ring finger and CHY zinc finger domain containing 1
GENE ID	25898
ORIGIN SPECIES	Human
ACCESSIONS	BC047393
VECTOR	pDONR223
NCBI SYMBOL	RCHY1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgacggcccgggaagatggcgccagcgggtcaagagcgaggtcagcggggctgcgagcactatgacagaggatgtctctaaaggccacctgtgtgacaagctttatactgcccgttgtgtcatgataacaatgaagatcatcaactagatcgctttaagtgaaggaagtgcagtgacataaactgtgaaaaaattcaacatgcccaacagactgtgagaatgtagcacattgttggagaataattatgcgatataatgccattgtttgacaaagataagaagcagatcactgtgaaaactgtggaaattgttaggattgtgcaagggaagattttccattgttgaaatgtaactatgcctagctatgaattcctcaaggagacacaagtgattgaaaatgtgtcccaacagaattgtccaataatgttggaggacattcacacatcccgtgtgtgtcctatgtcttccatgtggacatctttacatagaacgtgttatgaagaaatgttgaagaaggctacagatgtccattatgtatgcactctgtttagatgaccaggtattggagacagctggatgatgaagtagcacagactcctatgccatcagaatcagaacatgactgtggatattctctcaatgactgtaatggacgatccactgttcagttcatatattagcatgaaatgaagattgtgaatcctataactgtctaa gctggaggacgtagaatttcactggatcagcaatga