

TBC1D26-TBC1 domain family, member 26 Gene

Reference: PTXBC047400

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TBC1 domain family member 26		
For research use only.			

GENE DESCRIPTION	TBC1 domain family, member 26
GENE ID	353149
ORIGIN SPECIES	Human
ACCESSIONS	BC047400
VECTOR	pDONR223
NCBI SYMBOL	TBC1D26
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagatggatgggacccgtataacctgctgccaggggcaagcgaatcatcattactaagtatgagcaggacacccagctggggcagcagtggaacttggggcatgagcaggttgatgtcagaaaatacaccaataacctcgggattgtcatgagatggagctgccccgcgt
cagtgccctggagggtgaagcaagacgcaaggaaagtaaacgtaccaacaagtggcaaaagatgctgagactggacaaaatataggagcaccaagaagctgtctcaagagtatacaaaagtcattccctggcggfacggggccggcgctgacttttctagatattgacagaa
tcaagtccagaaccagcgaatataaggtcatgaaggagaaggcgaagaggtcctccagaatcatccactgcatccagctagatgtcagccacacctgcagaaacacatgatgttcatacaagattcggagtcaggcagcaggaattatgtgacatcctcgtggcctattctgcata
aacctgaggtgggtaccacagggaactgagccgcatcactgcatcctcctctgtgtctgagagggaagatgtttctggcgcttaccagttgctcgtggtgagaggcactccctgtgtattctacagcccaatactgcctggctcagaagctcctatgcaccaggagcag
gtgtcgcacaagtcctcccaagatcatga