

UBXN10-UBX domain protein 10 Gene

Reference: PTXBC047428

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UBXD3; UBX domain-containing protein 10; UBX domain containing 3; UBX domain-containing protein 3; UBX domain protein 10		

For research use only.

GENE DESCRIPTION	UBX domain protein 10
GENE ID	127733
ORIGIN SPECIES	Human
ACCESSIONS	BC047428
VECTOR	pDONR223
NCBI SYMBOL	UBXN10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccacagaagccctgtgaatatagcaccacctgagtgtagcactgtgtgtcagcacagcagtgacagcctcatttggcagccaaactcactaaatgacatgataaggcccaagtcgccaaggacggacaagaccgagctgcagaaatccaggcggtggaggtgtgcgc
tcacatataccatctccgctccagccattccctatgagttgccaagcagccaaaaccaggagcctgtgcacccaaatctccaaaccaggagcttctgatgagatccctgagctgcagcagcaagtaccactggggcttctctctcacaataagatccagtccttccatcaaca
gaaagaacctggaggaggagctgtggaaccgttgccaaaaggccagctcactgcaactgagcagatcctgggtctcttaccagacgagacgggcacatgaagacaagtgaaagaattccagagctcgagctgtgcccgtggagaggaaatcatctccgaaccaagaac
agggtcttccagggtggaaatctggaggaaacctcgacacgaaccaaaggtgtgtgtgtatgacacacagccaaagggtgtgtacgcatctccgccaacagatgattgcaaacattgttctgtgcccgaacagaaaacaaacctctaccgacactgcagcatt
gaaacaatggaggtgccaggaggcgattttctgacctcaccaaatctctgcaagatgcagaatccccacaaagctgtgtgtggcactcactggaagatggggaagggtggcctga