

TMEM17-transmembrane protein 17 Gene

Reference: PTXBC047439

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 17		
For research use only.			

GENE DESCRIPTION	transmembrane protein 17
GENE ID	200728
ORIGIN SPECIES	Human
ACCESSIONS	BC047439
VECTOR	pDONR223
NCBI SYMBOL	TMEM17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagctgccagatccggtgcgccagcgggtgggaaacttcagccggcggtgttcagtgattccaatcggaccagtcagagtcgaatgaggggtccggaatgaaatggtctcagttggcactgcagatgtcactttatttatactactatttccactgtgggtgggtgagcagca
ttatgatgttcacatgaagtattcaatctgtcctgactactacaaatcatgtgatcactgttatcatcctaataacctaattgaagccatccggtgttatctgggtacgtgggttaacctacaggagaagggttcctgagttggctggcttttggcttttgagccttctatgcagttacctaattctttt
ctgtctttaaagaaggcctaacaatctgcccttgaaaaagcgatacatatcatcttcactctctccttgccttccaagttgttcagcatttcttaccttaaggaaatgggttaacagttggcagttcgtttccactccaagacttgaccggctctctgcaaacagaggagacatgagaagg
atgaggtcatgtatagaagagatctga