

ANGEL2-angel homolog 2 (Drosophila) Gene

Reference: PTXBC047469

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Ccr4d; KIAA0759L; protein angel homolog 2; angel homolog 2		

For research use only.

GENE DESCRIPTION	angel homolog 2 (Drosophila)
GENE ID	90806
ORIGIN SPECIES	Human
ACCESSIONS	BC047469
VECTOR	pDONR223
NCBI SYMBOL	ANGEL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcctataatactttcacagaatttactggaagataactctacattatagacattgccggcggccagttattacactggagtttttaggtttcccaattattctgaagaaattaacattttgatgcagacgtactttgttgaagaagttcaagaagatcattatggagcagagatcaggccaa
gtttggaatcactgggttatcactgtgaatataagatcgagcaggaaggaaacctgatggctgtgctattttcttcaaacattccaaattttactctgtcagtgaaacccagtggaattctccgccctgataattctctgttggacagagacaattgttgatttttactctacagccaaaatt
ccatatgtcgtcctgccctgcaatctgcgtagcaaatagcactctgtgtataatccaaggcgagggtgatattaaagctgacgcaattggcaatgctactggcagagattccagtggtgccaccagaaagatggcagctctgccctattgttatgtgtgtgactttaattctgttctctcca
ctatatagttcataaaggaaagaaattgaattatgaaggacttccataggaaggatctgtggccaggaaacagctctcagggggacaaagaattttatctattccaatttggccccaacctaggtatctcacagaactgtgtgtatgaggtacagcaggtacaaaagtagaaagacaga
cagtgatctgacacaaacacagctgaagcaaacagaggtcctagtgcagctgaaaaattgtcttcaaatctacagcaccatttcagttgtcatctgtttattaa