

NXPH1-neurexophilin 1 Gene

Reference: PTXBC047505

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NPH1; Nbla00697; neurexophilin-1; neurexophilin 1		
For research use only.			

GENE DESCRIPTION	neurexophilin 1
GENE ID	30010
ORIGIN SPECIES	Human
ACCESSIONS	BC047505
VECTOR	pDONR223
NCBI SYMBOL	NXPH1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggctgcgtggtacgtctttctctcgcagccaccgtctctgtgacatgtgccattaacgaacggtggaaagtgcagaactctgaaatcaggagcagcaaatccacactaaagcacatattggacagaaagcagcaaaagactgtctatcagccgactcgtcacag
actttctgtggcaaaagagaatgatacagatttggacctgagatagacacccagaaccttattctgagcaagacctctgggactggctgaggaaactccacagacctcaagagcctcggccagggccaagagaaggccattgttaaacgggcaagttaagaaaatttggatggg
gctgatttattccaacatcaaaacagtgaagctgaacctgttgataactgggaaattgtatcatgtgcaatgggacatttagtatttattcaggcataattcaactggcaagggaatgtatctgtcagctgtgtacccctcaaaaaatcgtggaatttgacttggcacaacaaccgtgatt
gatgccaaagattccaagcttttaattgtcgactgaatgaaaaggttgacaaggctaccaagaacacactctgcaactatgacccctcaaaaacctgttaccaggagcaaacccaagtcattgtatcctgctctgtccaaagcccttaagggtgatctgtattacatttctttatagfaca
gattataaactgttacagaaggtgtgcctgactacaactaccacagtgacacaccttacttccctcgggatga