

DMRTC1B-DMRT-like family C1B Gene

Reference: PTXBC047596

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	doublesex- and mab-3-related transcription factor C1; doublesex-mab-3; DMRT like family C1		

For research use only.

GENE DESCRIPTION	DMRT-like family C1B
GENE ID	728656
ORIGIN SPECIES	Human
ACCESSIONS	BC047596
VECTOR	pDONR223
NCBI SYMBOL	DMRTC1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgcccctccaaagctcccatcgtgtcaggaatttgaccatcagagcaggagccctcactgggaaggagaaacaacatgtgcagcccgagaccacatcttcacagccccgaggaggtcccaagtctctgaccaggtttgtttctgcccaactcagagtgccagcggaa
gctggaggccgctgaggtctgtgactctgagaaactctgccaggccctcctgactccatctcctgcaccagccttgcaaccaccagctcctgctggagataaaggattccagctccagccctctctcgcgccagggccagccagctccatctcgtgcctattggacatctgg
gatgcattcctcttgagctag