

ZNF767-zinc finger family member 767 Gene

Reference: PTXBC047675

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger family member 767		
For research use only.			

GENE DESCRIPTION	zinc finger family member 767
GENE ID	79970
ORIGIN SPECIES	Human
ACCESSIONS	BC047675
VECTOR	pDONR223
NCBI SYMBOL	ZNF767
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggaggcggcgcagctccgatttctccgtggacgatggcagccacgattcagccatggagaggaagattgaatcgaggctgctcacctgtttccctagaagggtcaaaccgggatggccgagaagaagctggctgattgcgagaagacagccgtggagttcgggaaccag
ctggaggcgaagtggccgtgctggggacccctgctgcagagtagtggcctgctgcagagggcgtggagaacgtggagaacctgctgcgaacaggaaacttctggatcctgcggctgccccaggcagcaagggggagtccttaagggtgcctatgacattgatgatgtgctgtg
tag