

IL7-interleukin 7 Gene

Reference: PTXBC047698

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	interleukin-7; interleukin 7		

For research use only.

GENE DESCRIPTION	interleukin 7
GENE ID	3574
ORIGIN SPECIES	Human
ACCESSIONS	BC047698
VECTOR	pDONR223
NCBI SYMBOL	IL7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttccatgttttttaggtatattcttggacttcctccctgacctgtgtctgtgccagtagcatcatctgattgtgatattgaaggtaaagatggcaacaatatgagagtggtctaatggcagcatgatcaattattggacagcatgaaagaaattgtagcaattgcctgaataatgaattta
acttttttaaaagacatctgtgatgtaataaggaggtatgttttattccgtgctcgcgaagttgaggcaatttcttaaaatgaatgacactgggtattttgatctccactattaaaagttcagaaggcacaacaatactgtgaactgcactggccagggttaaaggagaaaaccagctgc
cctgggtgaagcccaaccaacaagagtttggagaaaataaactttaaaggaacagaaaaaactgaatgactgtgttcttaagagactattacaagataaaaactgttggataaaaattttgatgggcactaagaacactga