



Reference: PTXBC047719

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane channel-like protein 7; transmembrane channel like 7		
For research use only.			

GENE DESCRIPTION	transmembrane channel-like 7
GENE ID	79905
ORIGIN SPECIES	Human
ACCESSIONS	BC047719
VECTOR	pDONR223
NCBI SYMBOL	TMC7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr