



Reference: PTXBC047720

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DFNB77; LH2D1; lipoxygenase homology domain-containing protein 1; lipoxygenase homology domains 1		

For research use only.

GENE DESCRIPTION	lipoxygenase homology domains 1
GENE ID	125336
ORIGIN SPECIES	Human
ACCESSIONS	BC047720
VECTOR	pDONR223
NCBI SYMBOL	LOXHD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr