

C2orf64-chromosome 2 open reading frame 64 Gene

Reference: PTXBC047722

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein C2orf64; C2orf64; 6330578E17Rik; CEMCOX3; Pet191; cytochrome c oxidase assembly factor 5		

For research use only.

GENE DESCRIPTION	chromosome 2 open reading frame 64
GENE ID	493753
ORIGIN SPECIES	Human
ACCESSIONS	BC047722
VECTOR	pDONR223
NCBI SYMBOL	C2orf64
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctaagtattatgaggacaagccgagggcgccgctgcgcggcgctgaaggaggacgtggcgcgctgtctgctgcagtcggactgtgtgtccagggaaggaaaatcacctcggcagtggttgaagggaaggatactgcaactctttgaagtacgcatttttgagtgtaaaagatcag
tgttggataacagggcaagattcagaggagaaaaaggatattga