



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

OAF-OAF homolog (Drosophila) Gene

Reference: PTXBC047726

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

OAF homolog; NS5ATP13TP2; out at first protein homolog; HCV NS5A-transactivated protein 13 target protein 2; out at first homolog

For research use only.

GENE DESCRIPTION	OAF homolog (Drosophila)
GENE ID	220323
ORIGIN SPECIES	Human
ACCESSIONS	BC047726
VECTOR	pDONR223
NCBI SYMBOL	OAF
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr