

ABCA8-ATP-binding cassette, sub-family A (ABC1), member 8 Gene

Reference: PTXBC047765

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP-binding cassette sub-family A member 8; ATP-binding cassette, sub-family A (ABC1), member 8; ATP binding cassette subfamily A member 8		

For research use only.

GENE DESCRIPTION	ATP-binding cassette, sub-family A (ABC1), member 8
GENE ID	10351
ORIGIN SPECIES	Human
ACCESSIONS	BC047765
VECTOR	pDONR223
NCBI SYMBOL	ABCA8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggagagaaagatcagtggtgtcaacaactggccttattatgaagaactttcttaaaaaatggagaatgaaagagagatcctaatggatggctgaattcattgctcctactactttgtttgtatatatctcatagtcacagtaaatgattttctcactgcttaccatggacgtg
ggacgggtagatcatttaataatccagattttctgtgtatcacacctgtcaccaacacgacccaacagataatgaataaagtagcctctactcctctcgtggcaggtaaagaggctctgggactgccagatgaggaaagtattaaagaattcacagcaaatatcctgaagaaatagtaag
agtcacctttactaatacactacatcatttgaagtcttctagacatggatgccagcaagaaggagcacaaggaccatacagctcattgttatgaaacaaatgaagatgtttactgtgaagtttcagttttggaaggagggtttgtggctctcactgccaattatgctgctattata
gaaatcacacaatacactcagtgatggaggagctgatgtcagttactggaaaaaatatgaagatgcattcctcattggtcaatcaggagtataactgattgtacctttttcctgcattatttcatttctcatttactatgcattctgttaatgtcacaagagaggaaaaggatgaaggc
cttgatgacaatgatgggtcttcgggattcagcgttctgtga