

C11orf46-chromosome 11 open reading frame 46 Gene

Reference: PTXBC047775

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C11orf46; ARF7EP; dJ299F11.1; ARL14 effector protein; ADP-ribosylation factor-like 14 effector protein; ARF7 effector protein; ADP ribosylation factor like GTPase 14 effector protein		

For research use only.

GENE DESCRIPTION	chromosome 11 open reading frame 46
GENE ID	120534
ORIGIN SPECIES	Human
ACCESSIONS	BC047775
VECTOR	pDONR223
NCBI SYMBOL	C11orf46
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggatccatgttcagttggagtcagcttcgtactacaaatgagtgccataaaacctactactctgcacacaggttttaagactttgcaagaattgtcatcaaatgatatgcttttactcaactagaactggaatgacactttctgggaacaatacaattgctttcatatgtaaaatttaca
ttgacagatttgaggatttacagaagtcattgtgtgacccatttaacatacacacaagaattagccaaaaaatttgcattgtaattgacttagatgatgccacttttctgagtgctaaatttggaagacagctgtacctggttggaagctttgtccaaatgcacacagataatcaatggaagtgtgg
atgttgatactgaagaccgccagaaaaggaaacctgagtcagatggagaactgtctaaagctttgaggtcattacaatttacgaatccaggaaggcaaaactgaatttgctccagaactggtaaaagagaaaaaagagcttacaaaaatgcaaccgctggttcagacagacaagtgat
accagcaagagtaaggtctatgatagccagggtcctcctgatttttagtgggatggacctctgtgactgcctggatgaagactgcttaggatgtttctatgctgtcctgctgtggttctaccaagtgtggagctgaatccgctgtgaccgcaagtggctgtatgagcaaatgaaattgaagg
aggagaataattcataataaactgctggataa